METER

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USP 922

TYPES OF WATER ACTIVITY INSTRUMENTS

Dew Point Temperature—Chilled Mirror

At the dew point temperature, the relative humidity will be 100%, the evaporation and condensation rates will be equal, and there will not be any net condensation. When the temperature of a surface drops below the dew point temperature, water will condense on it. The dew point temperature can be determined by cooling a mirror in the sample chamber until condensation begins to occur and recording the mirror temperature at that point. The condensation on the chilled mirror is detected by an optical sensor. The difference between the dew point temperature and the temperature of the sample chamber is used to determine the relative humidity. The mirror must be kept clean and protected from contamination by dust or residue that may condense on the chilled mirror. Therefore, measurements on microleak patches may increase the likelihood of repeat calibration issues. The measurement of dew point temperature is a primary method of determining the vapor pressure of the headspace.

Electronic Hygrometers—Resistive (Electrolytic) Hygrometer

In this instrument, the relative humidity in the sample chamber is measured by a hygroscopic resistor. The resistor is a liquid electrolyte that absorbs or loses water to the air in the sample chamber, and the resulting resistance is proportional to the relative humidity. These resistors provide good precision and typically function over the a_w range of 0.04–1.00. Soluble substances, which are soluble in the liquid electrolyte, may affect the resistance and interfere with the relative humidity measurement. Chemical filters may be used to absorb volatiles and prevent interference. The response of these instruments is not necessarily linear, so sensors should be calibrated using multiple reference solutions over the entire operating range. Alternatively, a two-point linear calibration over a narrower range that brackets the sample may be used.

Electronic Hygrometers—Capacitive Hygrometer

In this instrument, the relative humidity in the sample chamber is measured by a hygroscopic capacitor. The capacitor consists of a hygroscopic polymer held between conductive plates. The polymer absorbs or loses water to the air in the sample chamber, and the resulting capacitance is proportional to the relative humidity. These sensors may become saturated at high a_w and may exhibit hysteresis after measurement of high a_w samples. Capacitive hygrometers are typically less affected by water vapors than resistive or dew point sensors; however, care should still be used when performing measurements in the presence of other volatile compounds. The response of these instruments is not necessarily linear, so sensors should be calibrated using multiple reference solutions over the entire operating range. Alternatively, a two-point linear calibration over a narrower range that brackets the sample may be used.

Optical Hygrometer—Tunable Diode Laser

Water molecules in the gas phase absorb a product about 1/3 and 1/34 radiation. Strong and broad absorption bands exist that can be used to measure the vapor pressure of air, but many of these are near absorption bands of other volatiles that may interfere. The most specialized optical hygrometers are not used for water activity measurements. The absorption bands are composed of individual transitions and unique not only to a specific molecule but to a particular isotope of water. A diode laser can produce light at a specific wavelength, with a narrow bandwidth of ± 1 nm. The wavelength of a tunable diode laser (TDL) can be tuned to a specific, isolated absorption line in the water vapor spectrum by controlling the laser temperature and can be scanned across that line by modulating the laser current (and, therefore, the laser's temperature). This modulation provides both the absorption and unperturbed laser intensity in the headspace, from which the vapor pressure of the headspace can be calculated. Because the TDL is capable of measuring vapor pressure in the presence of even high concentrations of other volatiles, it is the only method that can correctly measure water activity in these situations. Electronic and dew point hygrometers are sometimes suitable in the presence of low concentrations of volatiles, but readings are always adversely affected by the volatiles. The TDL also can be used to measure water activity noninvasively in sealed and sealed containers.

Near-Infrared Spectroscopy

Near-infrared spectroscopy (NIR) data has traditionally been used to make the a_w measurements with the water content. The use of NIR to measure a_w is based on the generation of a predictive model for a particular official article (product or raw material) to correlate NIR absorbance spectra to the a_w measured with one of the previously mentioned direct methods. This predictive model may be used to rapidly predict the a_w of additional samples. When additional predictive models are available, NIR can be used for simultaneous determination of other compounds such as the content of active ingredients, additives, or impurities. The NIR method is best suited for rapid analysis of raw material, in-process materials, and end products.

Other Methods

Other methods for measurement of a_w may be useful but are not currently employed for routine measurements of pharmaceutical products, dietary supplements, and raw materials. These methods include gravimetry, lowing point depression, thermogravimetric analysis, and aqueous methods. The hygrometers may be used for approximate a_w measurements but not for routine work due to very high precision times and low accuracy. The freezing-point and thermocouple

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